


# ORDERING INFORMATION

• Based on Delivery Device : Distal  Proximal

- (NNN) Non covered
- (CCC) Fully covered
- (NCN) Partially covered

XXX-00-000-000 XXX 00 000-X000  
 ① ② ③ ④ ① ② ③ ④  
 ① Product Name ③ Stent Length  
 ② Stent Diameter ④ Delivery Device Length



## ESOPHAGUS




### HANAROSTENT® Esophagus Upper (CCC)



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
ECP18060-X070 ⋮ ECP18170-X070	24-18-26	32~142*	60~170*	700	6/18
ECP20060-X070 ⋮ ECP20170-X070	26-20-28	32~142*	60~170*	700	6/18

\*Increment of 10mm


### CHOOSTENT® Esophagus Upper (CCC)



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
NESC-18-060-070 ⋮ NESC-18-170-070	24-18-26	32~142*	60~170*	700	6/18

\*Increment of 10mm

### HANAROSTENT® Esophagus Asymmetric (CCC)



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
ECN18060-X070 ⋮ ECN18170-X070	24-18-24	30~140*	60~170*	700	6/18
ECN20060-X070 ⋮ ECN20170-X070	26-20-26	30~140*	60~170*	700	6/18
ECN24060-Z070 ⋮ ECN24170-Z070	30-24-30	30~140*	60~170*	700	8/24

\*Increment of 10mm

### CHOOSTENT® Esophagus Asymmetric (CCC)



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
NESA-18-060-070 ⋮ NESA-18-170-070	24-18-24	30~140*	60~170*	700	6/18
NESA-24-060-070 ⋮ NESA-24-170-070	30-24-30	30~140*	60~170*	700	8/24

\*Increment of 10mm

### HANAROSTENT® Esophagus (CCC)



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
HES-18-060-070 ⋮ HES-18-170-070	24-18-24	20~130*	60~170*	700	6/18
HES-20-060-070 ⋮ HES-20-170-070	26-20-26	20~130*	60~170*	700	6/18
HES-22-060-070 ⋮ HES-22-170-070	28-22-28	20~130*	60~170*	700	8/24

\*Increment of 10mm


### CHOOSTENT® Esophagus (CCC)



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
NES-18-060-070 ⋮ NES-18-170-070	24-18-24	20~130*	60~170*	700	6/18
NES-20-060-070 ⋮ NES-20-170-070	26-20-26	20~130*	60~170*	700	6/18
NES-22-060-070 ⋮ NES-22-170-070	28-22-28	20~130*	60~170*	700	8/24

\*Increment of 10mm


### HANAROSTENT® Esophagus Valve (CCC)



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
HESV-18-060-070 ⋮ HESV-18-170-070	24-18-24	20~130*	60~170*	700	8/24
HESV-22-060-070 ⋮ HESV-22-170-070	28-22-28	20~130*	60~170*	700	8/24

\*Increment of 10mm

### CHOOSTENT® Esophagus Valve (CCC)



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
HEV-18-060-070 ⋮ HEV-18-170-070	24-18-24	20~130*	60~170*	700	8/24
HEV-22-060-070 ⋮ HEV-22-170-070	28-22-28	20~130*	60~170*	700	8/24

\*Increment of 10mm

**CHOOSTENT®**  
**Esophagus Upper ST (CCC)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
NESCT-18-060-070 ⋮ NESCT-18-170-070	24-18-24	32~142*	60~170*	700	8/24

\*Increment of 10mm

**CHOOSTENT®**  
**Esophagus ST (CCN)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
EPG18060-X070 ⋮ EPG18170-X070	24-18-24	20~130*	60~170*	700	6/18
EPG22060-Z070 ⋮ EPG22170-Z070	28-22-28	20~130*	60~170*	700	8/24

\*Increment of 10mm

**CHOOSTENT®**  
**Esophagus Valve ST (CCN)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
HEVT-18-060-070 ⋮ HEVT-18-170-070	24-18-24	20~130*	60~170*	700	8/24
HEVT-22-060-070 ⋮ HEVT-22-170-070	28-22-28	20~130*	60~170*	700	8/24

\*Increment of 10mm

**HANAROSTENT®**  
**Esophagus BS (CCC)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
ECV20060-X070 ⋮ ECV20170-X070	26-20-26	20~130*	60~170*	700	6/18
ECV22060-Z070 ⋮ ECV22170-Z070	28-22-28	20~130*	60~170*	700	8/24

\*Increment of 10mm


**HANAROSTENT®**  
**Esophagus BS (CCC)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
NES-28-060-070 ⋮ NES-28-150-070	34-28-34	20~110*	60~150*	700	8/24

\*Increment of 10mm


**HANAROSTENT®**  
**Esophagus Valve BS (CCN)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
EPV-18-060-070 ⋮ EPV-18-170-070	24-18-24	20~130*	60~170*	700	8/24
EPV-22-060-070 ⋮ EPV-22-170-070	28-22-28	20~130*	60~170*	700	8/24

\*Increment of 10mm

**HANAROSTENT®**  
**Esophagus Valve BS (CCN)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
ECW18060-Z070 ⋮ ECW18170-Z070	24-18-24	20~130*	60~170*	700	8/24
ECW22060-Z070 ⋮ ECW22170-Z070	28-22-28	20~130*	60~170*	700	8/24

\*Increment of 10mm

**HANAROSTENT®**  
**Esophagus BS (NCN)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
EPW-20-060-070 ⋮ EPW-20-170-070	26-20-26	20~130*	60~170*	700	6/18

\*Increment of 10mm

**HANAROSTENT®**  
**Esophagus (NCN)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
EPB18060-X070 ⋮ EPB18170-X070	24-18-24	20~130*	60~170*	700	6/18
EPB22060-Z070 ⋮ EPB22170-Z070	28-22-28	20~130*	60~170*	700	8/24

\*Increment of 10mm


**CHOOSTENT®**  
**Esophagus (NCN)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
EPE18060-X070 ⋮ EPE18170-X070	24-18-24	20~130*	60~170*	700	6/18
EPE20060-X070 ⋮ EPE20170-X070	26-20-26	20~130*	60~170*	700	6/18
EPE22060-Z070 ⋮ EPE22170-Z070	28-22-28	20~130*	60~170*	700	8/24

\*Increment of 10mm


**HANAROSTENT®**  
**Esophagus Valve (CCN)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
EPC-18-060-070 ⋮ EPC-18-170-070	24-18-24	20~130*	60~170*	700	8/24
EPC-22-060-070 ⋮ EPC-22-170-070	28-22-28	20~130*	60~170*	700	8/24

\*Increment of 10mm

**CHOOSTENT®**  
**Esophagus Valve (CCN)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
EPF18060-Z070 ⋮ EPF18170-Z070	24-18-24	20~130*	60~170*	700	8/24
EPF22060-Z070 ⋮ EPF22170-Z070	28-22-28	20~130*	60~170*	700	8/24

\*Increment of 10mm

**HANAROSTENT®**  
**Esophagus TTS (CCC)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
ECBA-18-060-180 ⋮ ECBA-18-150-180	26-18-26	30-120*	60-150*	1800	3.5/10.5
ECBA-18-060-230 ⋮ ECBA-18-150-230	26-18-26	30-120*	60-150*	2300	3.5/10.5

\*Increment of 10mm

**HANAROSTENT®**  
**Esophagus TTS (NCN)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
EPBA-18-060-180 ⋮ EPBA-18-150-180	26-18-26	30-120*	60-150*	1800	3.5/10.5
EPBA-20-060-180 ⋮ EPBA-20-150-180	26-20-26	30-120*	60-150*	1800	3.5/10.5
EPBA-18-060-230 ⋮ EPBA-18-150-230	26-18-26	30-120*	60-150*	2300	3.5/10.5
EPBA-20-060-230 ⋮ EPBA-20-150-230	26-20-26	30-120*	60-150*	2300	3.5/10.5

\*Increment of 10mm

**HANAROSTENT®**  
**Esophagus Benign BS (CCC)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
EBN18060-Z070 ⋮ EBN18120-Z070	24-18-24	20-80*	60-120*	700	8/24
EBN20060-Z070 ⋮ EBN20120-Z070	26-20-26	20-80*	60-120*	700	8/24
EBN22060-Z070 ⋮ EBN22120-Z070	28-22-28	20-80*	60-120*	700	8/24

\*Increment of 20mm


**HANAROSTENT®**  
**Esophagus Benign BS ST (CCC)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
EBNT18060-Z070 ⋮ EBNT18120-Z070	24-18-24	20-80*	60-120*	700	8/24
EBNT20060-Z070 ⋮ EBNT20120-Z070	26-20-26	20-80*	60-120*	700	8/24
EBNT22060-Z070 ⋮ EBNT22120-Z070	28-22-28	20-80*	60-120*	700	8/24

\*Increment of 20mm

**HANAROSTENT®**  
**Esophagus Skidproof (CCC)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
ECD-20-080-070	26-20-26	40	80	700	6/18
ECD-20-100-070		60	100	700	6/18
ECD-20-120-070		80	120	700	6/18
ECD-20-140-070		100	140	700	6/18
ECD-20-160-070		120	160	700	6/18

**HANAROSTENT®**  
**Esophagus Lasso (NNN)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
END18060-C070 ⋮ END18170-C070	24-18-24	20-130*	60-170*	700	4/12
END20060-C070 ⋮ END20170-C070	26-20-26	20-130*	60-170*	700	4/12

\*Increment of 10mm

**HANAROSTENT®**  
**Esophagus Bariatric Surgery (CCC)**




Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
ECBB-28-180-090	32-28-28	145	180	900	8/24
ECBB-28-210-090		175	210	900	8/24
ECBB-28-240-090		205	240	900	8/24

**HANAROSTENT®**  
**Esophagus Bariatric Surgery ST (CCC)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
ECBBT-28-180-090	32-28-28	145	180	900	8/24
ECBBT-28-210-090		175	210	900	8/24
ECBBT-28-240-090		205	240	900	8/24

**HANAROSTENT®**  
**Gastro-Seal™ Esophagus (CCC)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
ECBS-26-180-090	36-26-34	130	180	900	8/24
ECBS-26-210-090		160	210	900	8/24
ECBS-26-240-090		190	240	900	8/24

**HANAROSTENT®**  
**Gastro-Seal™ Esophagus ST (CCC)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
ECBST-26-180-090	36-26-34	130	180	900	8/24
ECBST-26-210-090		160	210	900	8/24
ECBST-26-240-090		190	240	900	8/24

# BILIARY

## HANAROSTENT®

### Biliary (NNN)



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (Fr)
SHS-08-040-180 SHS-08-120-180	11.5-8-11.5	26~106*	40~120*	1800	7
SHS-10-040-180 SHS-10-120-180	13.5-10-13.5	26~106*	40~120*	1800	7
SHS-08-040-060 SHS-08-120-060	11.5-8-11.5	26~106*	40~120*	600	7
SHS-10-040-060 SHS-10-120-060	13.5-10-13.5	26~106*	40~120*	600	7

\*Increment of 10mm

## HANAROSTENT®

### FAST™ Biliary (NNN)



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (Fr)
SHS-08-040-S180 SHS-08-100-S180	11.5-8-11.5	26~86*	40~100*	1800	8.5
SHS-10-040-S180 SHS-10-100-S180	13.5-10-13.5	26~86*	40~100*	1800	8.5

\*Increment of 10mm

## HANAROSTENT®

### Biliary (NNN)




Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (Fr)
BNA08040-G180 BNA08120-G180	11.5-8-11.5	26~106*	40~120*	1800	8
BNA10040-G180 BNA10120-G180	13.5-10-13.5	26~106*	40~120*	1800	8
BNA08040-H060 BNA08120-H060	11.5-8-11.5	26~106*	40~120*	600	8
BNA10040-H060 BNA10120-H060	13.5-10-13.5	26~106*	40~120*	600	8

\*Increment of 10mm

## HANAROSTENT®

### Biliary (NNN)




Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (Fr)
BNL-08-040-180 BNL-08-100-180	11.5-8-11.5	26~86*	40~100*	1800	7
BNL-10-040-180 BNL-10-100-180	13.5-10-13.5	26~86*	40~100*	1800	7
BNL-08-040-060 BNL-08-100-060	11.5-8-11.5	26~86*	40~100*	600	7
BNL-10-040-060 BNL-10-100-060	13.5-10-13.5	26~86*	40~100*	600	7

\*Increment of 10mm

## HANAROSTENT®

### Biliary Hilar (NNN)



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (Fr)
BNY-08-060-180 BNY-08-120-180	11.5-8-11.5	46~106*	60~120*	1800	7
BNY-10-060-180 BNY-10-120-180	13.5-10-13.5	46~106*	60~120*	1800	7
BNY-08-060-060 BNY-08-120-060	11.5-8-11.5	46~106*	60~120*	600	7
BNY-10-060-060 BNY-10-120-060	13.5-10-13.5	46~106*	60~120*	600	7

\*Increment of 20mm

## HANAROSTENT®

### Benefit™ Biliary (NNN)




Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (Fr)
BNHS-06-040-180 BNHS-06-120-180	6	-	40~120*	1800	5.9
BNHS-08-040-180 BNHS-08-120-180	8	-	40~120*	1800	5.9
BNHS-06-040-060 BNHS-06-120-060	6	-	40~120*	600	5.9
BNHS-08-040-060 BNHS-08-120-060	8	-	40~120*	600	5.9

\*Increment of 10mm

## HANAROSTENT®

### Biliary Valve (NNN)




※ The stent is reversely loaded in percutaneous application

Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (Fr)
BNX-08-060-180 BNX-08-120-180	11.5-8	53~113*	60~120*	1800	10.2
BNX-10-060-180 BNX-10-120-180	13.5-10	53~113*	60~120*	1800	10.2
BNX-08-060-060 BNX-08-120-060	8-11.5	53~113*	60~120*	600	10.2
BNX-10-060-060 BNX-10-120-060	10-13.5	53~113*	60~120*	600	10.2

\*Increment of 10mm

## HANAROSTENT®

### Biliary (NC)




※ The stent is reversely loaded in percutaneous application

Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (Fr)
BPN-08-040-180 BPN-08-120-180	11.5-8-11.5	26~106*	40~120*	1800	8.5
BPN-10-040-180 BPN-10-120-180	13.5-10-13.5	26~106*	40~120*	1800	8.5
BPN-08-040-060 BPN-08-120-060	11.5-8-11.5	26~106*	40~120*	600	8.5
BPN-10-040-060 BPN-10-120-060	13.5-10-13.5	26~106*	40~120*	600	8.5

\*Increment of 10mm


**HANAROSTENT®**  
**Biliary (NCN)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (Fr)
BPA08040-D180 BPA08120-D180	11.5-8-11.5	26~106*	40~120*	1800	8.5
BPA10040-D180 BPA10120-D180	13.5-10-13.5	26~106*	40~120*	1800	8.5
BPA08040-R060 BPA08120-R060	11.5-8-11.5	26~106*	40~120*	600	8.5
BPA10040-R060 BPA10120-R060	13.5-10-13.5	26~106*	40~120*	600	8.5

\*Increment of 10mm

**HANAROSTENT®**  
**Biliary (NC)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (Fr)
BPD-10-060-180 BPD-10-100-180	10-20	54~94*	60~100*	1800	8.5

\*Increment of 10mm

**HANAROSTENT®**  
**Biliary Lasso (CCC)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (Fr)
SHCL-06-040-180 SHCL-06-120-180	8.3-6-8.3	30~110*	40~120*	1800	8.5
SHCL-08-040-180 SHCL-08-120-180	11.5-8-11.5	26~106*	40~120*	1800	8.5
SHCL-10-040-180 SHCL-10-120-180	13.5-10-13.5	26~106*	40~120*	1800	8.5
SHCL-06-040-060 SHCL-06-120-060	8.3-6-8.3	30~110*	40~120*	600	8.5
SHCL-08-040-060 SHCL-08-120-060	11.5-8-11.5	26~106*	40~120*	600	8.5
SHCL-10-040-060 SHCL-10-120-060	13.5-10-13.5	26~106*	40~120*	600	8.5

\*Increment of 10mm

**HANAROSTENT®**  
**Plumber™ (CCC)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (Fr)
BCF-10-020-180 BCF-10-040-180	22-10-22	13~33*	20~40*	1800	10.5
BCF-12-020-180 BCF-12-040-180	24-12-24	13~33*	20~40*	1800	10.5
BCF-14-020-180 BCF-14-040-180	26-14-26	12~32*	20~40*	1800	10.5
BCF-16-020-180 BCF-16-040-180	28-16-28	12~32*	20~40*	1800	10.5

\*Increment of 10mm

**HANAROSTENT®**  
**Biliary (CCC)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (Fr)
SHC-08-040-180 SHC-08-120-180	11.5-8-11.5	26~106*	40~120*	1800	8.5
SHC-10-040-180 SHC-10-120-180	13.5-10-13.5	26~106*	40~120*	1800	8.5
SHC-08-040-060 SHC-08-120-060	11.5-8-11.5	26~106*	40~120*	600	8.5
SHC-10-040-060 SHC-10-120-060	13.5-10-13.5	26~106*	40~120*	600	8.5

\*Increment of 10mm

**HANAROSTENT®**  
**Biliary Lasso (CCC)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (Fr)
BCB08040-E180 BCB08120-E180	11.5-8-11.5	30~110*	40~120*	1800	8
BCB10040-E180 BCB10120-E180	13.5-10-13.5	30~110*	40~120*	1800	8
BCB08040-Q060 BCB08120-Q060	11.5-8-11.5	30~110*	40~120*	600	8
BCB10040-Q060 BCB10120-Q060	13.5-10-13.5	30~110*	40~120*	600	8

\*Increment of 10mm

**HANAROSTENT®**  
**Biliary Lasso (CCC)**

※ Endoscopic application



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (Fr)
BCBC08040-B180 BCBC08120-B180	11.5-8-11.5	26~106*	40~120*	1800	10.2
BCBC10040-B180 BCBC10120-B180	13.5-10-13.5	26~106*	40~120*	1800	10.2

※ Percutaneous application



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (Fr)
BCBC08040-N060 BCBC08120-N060	11.5-8-11.5	26~106*	40~120*	600	10.2
BCBC10040-N060 BCBC10120-N060	13.5-10-13.5	26~106*	40~120*	600	10.2

\*Increment of 10mm

**HANAROSTENT®**  
Multi-Hole™ Biliary (CCC)



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (Fr)
BMCL-08-040-180 ⋮ BMCL-08-120-180	11.5-8-11.5	26~106*	40~120*	1800	8
BMCL-10-040-180 ⋮ BMCL-10-120-180	13.5-10-13.5	26~106*	40~120*	1800	8
BMCL-08-040-060 ⋮ BMCL-08-120-060	11.5-8-11.5	26~106*	40~120*	600	8
BMCL-10-040-060 ⋮ BMCL-10-120-060	13.5-10-13.5	26~106*	40~120*	600	8

\*Increment of 10mm

**HANAROSTENT®**  
Biliary Flap (CCC)




※ The stent is reversely loaded in percutaneous application

Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (Fr)
BCS08040-D180 ⋮ BCS08120-D180	8-11.5	25~105*	40~120*	1800	8.5
BCS10040-D180 ⋮ BCS10120-D180	10-13.5	25~105*	40~120*	1800	8.5
BCS08040-R060 ⋮ BCS08120-R060	11.5-8	25~105*	40~120*	600	8.5
BCS10040-R060 ⋮ BCS10120-R060	13.5-10	25~105*	40~120*	600	8.5

\*Increment of 10mm

**HANAROSTENT®**  
Biliary Flap Lasso (CCC)



※ Endoscopic application

Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (Fr)
BCT-08-040-180 ⋮ BCT-08-120-180	8-11.5	23~103*	40~120*	1800	9
BCT-10-040-180 ⋮ BCT-10-120-180	10-13.5	23~103*	40~120*	1800	9

※ Percutaneous application



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (Fr)
BCT-08-040-060 ⋮ BCT-08-120-060	8-11.5	23~103*	40~120*	600	9
BCT-10-040-060 ⋮ BCT-10-120-060	10-13.5	23~103*	40~120*	600	9

\*Increment of 10mm

**HANAROSTENT®**  
Biliary Flap (CCC)



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (Fr)
BCG-06-040-180 ⋮ BCG-06-120-180	6	20~100*	40~120*	1800	8
BCG-08-040-180 ⋮ BCG-08-120-180	8	20~100*	40~120*	1800	8.5
BCG-10-040-180 ⋮ BCG-10-120-180	10	20~100*	40~120*	1800	8.5
BCG-06-040-060 ⋮ BCG-06-120-060	6	20~100*	40~120*	600	8
BCG-08-040-060 ⋮ BCG-08-120-060	8	20~100*	40~120*	600	8.5
BCG-10-040-060 ⋮ BCG-10-120-060	10	20~100*	40~120*	600	8.5

\*Increment of 10mm

**HANAROSTENT®**  
FAST™ Biliary Flap (CCC)




Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (Fr)
BCG-08-040-S180 ⋮ BCG-08-080-S180	8	20~60*	40~80*	1800	8.5
BCG-10-040-S180 ⋮ BCG-10-080-S180	10	20~60*	40~80*	1800	8.5

\*Increment of 10mm

**CHOOSTENT®**  
**Colon/Rectum (CCC)**


Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
HCC-22-060-070 HCC-22-170-070	30-22-30	20~130*	60~170*	700	8/24
HCC-24-060-070 HCC-24-170-070	32-24-32	20~130*	60~170*	700	8/24

\*Increment of 10mm

**HANAROSTENT®**  
**Colon/Rectum (CCC)**


Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
CCB22060-Z070 CCB22170-Z070	30-22-30	20~130*	60~170*	700	8/24
CCB24060-Z070 CCB24170-Z070	32-24-32	20~130*	60~170*	700	8/24

\*Increment of 10mm

**HANAROSTENT®**  
**Colon/Rectum TTS (CCC)**


Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
CCBA-20-060-230 CCBA-20-140-230	26-20-26	24~104*	60~140*	2300	3.4/10.2

\*Increment of 10mm

**HANAROSTENT®**  
**Intestine (CCC)**


Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
HRC-16-040-230 HRC-16-090-230	26-16-26	24~74*	40~90*	2300	3.4/10.2
HRC-16-040-140 HRC-16-090-140	26-16-26	24~74*	40~90*	1400	3.4/10.2

\*Increment of 10mm

**HANAROSTENT®**  
**Colon/Rectum (NNN)**


Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
NCS-22-060-230 NCS-22-170-230	28-22-28	30~140*	60~170*	2300	3.4/10.2
NCS-24-060-230 NCS-24-170-230	30-24-30	30~140*	60~170*	2300	3.4/10.2
NCS-22-060-140 NCS-22-170-140	28-22-28	30~140*	60~170*	1400	3.4/10.2
NCS-24-060-140 NCS-24-170-140	30-24-30	30~140*	60~170*	1400	3.4/10.2
NCS-26-060-140 NCS-26-150-140	32-26-32	30~120*	60~150*	1400	4/12

\*Increment of 10mm

**CHOOSTENT®**  
**Colon/Rectum Asymmetric (CCC)**


Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
CCI22060-Z070 CCI22170-Z070	30-22-30	30~140*	60~170*	700	8/24
CCI24060-Z070 CCI24170-Z070	32-24-32	30~140*	60~170*	700	8/24

\*Increment of 10mm

**HANAROSTENT®**  
**Colon/Rectum Asymmetric (CCC)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
CCN22060-Z070 CCN22170-Z070	30-22-30	30~140*	60~170*	700	8/24
CCN24060-Z070 CCN24170-Z070	32-24-32	30~140*	60~170*	700	8/24

\*Increment of 10mm

**HANAROSTENT® TwynLay™**  
**Colon/Rectum (NCN)**


Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
TLC-20-060-230 TLC-20-150-230	26-20-26	30~120*	60~150*	2300	3.4/10.2
TLC-22-060-230 TLC-22-150-230	28-22-28	30~120*	60~150*	2300	3.4/10.2
TLC-20-060-140 TLC-20-150-140	26-20-26	30~120*	60~150*	1400	3.4/10.2
TLC-22-060-140 TLC-22-150-140	28-22-28	30~120*	60~150*	1400	3.4/10.2

\*Increment of 10mm

**HANAROSTENT®**  
**Intestine Asymmetric (NCN)**


Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
HRC-20-060-230 HRC-20-150-230	40-20-26	20~110*	60~150*	2300	3.4/10.2

\*Increment of 10mm

**HANAROSTENT®**  
**Colon/Rectum Lasso (NNN)**


Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
NCSL-22-060-230 NCSL-22-170-230	28-22-28	30~140*	60~170*	2300	3.4/10.2
NCSL-24-060-230 NCSL-24-170-230	30-24-30	30~140*	60~170*	2300	3.4/10.2
NCSL-22-060-140 NCSL-22-170-140	28-22-28	30~140*	60~170*	1400	3.4/10.2
NCSL-24-060-140 NCSL-24-170-140	30-24-30	30~140*	60~170*	1400	3.4/10.2

\*Increment of 10mm

**HANAROSTENT®****Colon/Rectum Lasso (NNN)**

Double layered Structure

Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
CNZ-22-060-230 CNZ-22-170-230	28-22-28	30~140*	60~170*	2300	3.4/10.2
CNZ-24-060-230 CNZ-24-170-230	30-24-30	30~140*	60~170*	2300	3.4/10.2
CNZ-22-060-140 CNZ-22-170-140	28-22-28	30~140*	60~170*	1400	3.4/10.2
CNZ-24-060-140 CNZ-24-170-140	30-24-30	30~140*	60~170*	1400	3.4/10.2
CNZ-26-060-140 CNZ-26-150-140	32-26-32	30~120*	60~150*	1400	4/12

\*Increment of 10mm

**HANAROSTENT® LowAx™****Colon/Rectum (NNN)**

Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
CNJF-18-060-230 CNJF-18-120-230	23-18	45~105*	60~120*	2300	3.4/10.2
CNJF-20-060-230 CNJF-20-120-230	25-20	45~105*	60~120*	2300	3.4/10.2
CNJF-22-060-230 CNJF-22-120-230	27-22	45~105*	60~120*	2300	3.4/10.2
CNJF-25-060-230 CNJF-25-120-230	30-25	45~105*	60~120*	2300	3.4/10.2
CNJF-18-060-140 CNJF-18-120-140	23-18	45~105*	60~120*	1400	3.4/10.2
CNJF-20-060-140 CNJF-20-120-140	25-20	45~105*	60~120*	1400	3.4/10.2
CNJF-22-060-140 CNJF-22-120-140	27-22	45~105*	60~120*	1400	3.4/10.2
CNJF-25-060-140 CNJF-25-120-140	30-25	45~105*	60~120*	1400	3.4/10.2

\*Increment of 10mm

**HANAROSTENT®**  
**Duodenum/Pylorus (NNN)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
NDS-18-060-230 NDS-18-170-230	24-18-24	30~140*	60~170*	2300	3.4/10.2
NDS-20-060-230 NDS-20-170-230	26-20-26	30~140*	60~170*	2300	3.4/10.2
NDS-18-060-120 NDS-18-170-120	24-18-24	30~140*	60~170*	1200	3.4/10.2
NDS-20-060-120 NDS-20-170-120	26-20-26	30~140*	60~170*	1200	3.4/10.2

\*Increment of 10mm

**HANAROSTENT®**  
**Duodenum/Pylorus Lasso (NNN)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
NDSL-18-060-230 NDSL-18-170-230	24-18-24	30~140*	60~170*	2300	3.4/10.2
NDSL-20-060-230 NDSL-20-170-230	26-20-26	30~140*	60~170*	2300	3.4/10.2
NDSL-18-060-120 NDSL-18-170-120	24-18-24	30~140*	60~170*	1200	3.4/10.2
NDSL-20-060-120 NDSL-20-170-120	26-20-26	30~140*	60~170*	1200	3.4/10.2

\*Increment of 10mm


**HANAROSTENT®**  
**Duodenum/Pylorus Lasso (NNN)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
DNZL-20-060-230 DNZL-20-170-230	26-20-26	30~140*	60~170*	2300	3.4/10.2
DNZL-22-060-230 DNZL-22-170-230	28-22-28	30~140*	60~170*	2300	3.4/10.2
DNZL-24-060-230 DNZL-24-170-230	30-24-30	30~140*	60~170*	2300	3.4/10.2
DNZL-20-060-120 DNZL-20-170-120	26-20-26	30~140*	60~170*	1200	3.4/10.2
DNZL-22-060-120 DNZL-22-170-120	28-22-28	30~140*	60~170*	1200	3.4/10.2
DNZL-24-060-120 DNZL-24-170-120	30-24-30	30~140*	60~170*	1200	3.4/10.2

\*Increment of 10mm


**HANAROSTENT®**  
**Duodenum/Pylorus (NCN)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
NDC-18-060-230 NDC-18-140-230	23-18-23	24~104*	60~140*	2300	3.4/10.2
NDC-20-060-230 NDC-20-140-230	26-20-26	24~104*	60~140*	2300	3.4/10.2
NDC-18-060-120 NDC-18-140-120	23-18-23	24~104*	60~140*	1200	3.4/10.2
NDC-20-060-120 NDC-20-140-120	26-20-26	24~104*	60~140*	1200	3.4/10.2

\*Increment of 10mm


**HANAROSTENT®**  
**Duodenum/Pylorus Lasso (NCN)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
NDCL-18-060-230 NDCL-18-140-230	23-18-23	24~104*	60~140*	2300	3.4/10.2
NDCL-20-060-230 NDCL-20-140-230	26-20-26	24~104*	60~140*	2300	3.4/10.2
NDCL-18-060-120 NDCL-18-140-120	23-18-23	24~104*	60~140*	1200	3.4/10.2
NDCL-20-060-120 NDCL-20-140-120	26-20-26	24~104*	60~140*	1200	3.4/10.2

\*Increment of 10mm

**HANAROSTENT®**  
**Duodenum/Pylorus Kim's Flare (NCN)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
DPDL-20-060-230 DPDL-20-150-230	22-20-40	32~122*	60~150*	2300	3.4/10.2
DPDL-20-060-120 DPDL-20-150-120	22-20-40	32~122*	60~150*	1200	3.4/10.2

\*Increment of 10mm

**HANAROSTENT® LowAx™**  
**Duodenum/Pylorus (NNN)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
DNJF-18-060-230 DNJF-18-120-230	18-23	45~105*	60~120*	2300	3.4/10.2
DNJF-20-060-230 DNJF-20-120-230	20-25	45~105*	60~120*	2300	3.4/10.2
DNJF-22-060-230 DNJF-22-120-230	22-27	45~105*	60~120*	2300	3.4/10.2
DNJF-18-060-120 DNJF-18-120-120	18-23	45~105*	60~120*	1200	3.4/10.2
DNJF-20-060-120 DNJF-20-120-120	20-25	45~105*	60~120*	1200	3.4/10.2
DNJF-22-060-120 DNJF-22-120-120	22-27	45~105*	60~120*	1200	3.4/10.2

\*Increment of 10mm

**HANAROSTENT® TwynLay™**  
**Duodenum/Pylorus (NCN)**



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length (mm)	Diameter (mm/Fr)
TLD-20-060-230 TLD-20-150-230	26-20-26	30~120*	60~150*	2300	3.4/10.2
TLD-22-060-230 TLD-22-150-230	28-22-28	30~120*	60~150*	2300	3.4/10.2
TLD-20-060-120 TLD-20-150-120	26-20-26	30~120*	60~150*	1200	3.4/10.2
TLD-22-060-120 TLD-22-150-120	28-22-28	30~120*	60~150*	1200	3.4/10.2

\*Increment of 10mm



## TRACHEOBRONCHIAL

### HANAROSTENT®

#### Trachea/Bronchium (CCC)



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(mm/Fr)
NTS-10-030-050 NTS-10-080-050	14-10-14	20~70*	30~80*	500 500	6/18 6/18
NTS-12-030-050 NTS-12-080-050	16-12-16	20~70*	30~80*	500 500	6/18 6/18
NTS-14-030-050 NTS-14-080-050	16.5-14-16.5	20~70*	30~80*	500 500	6/18 6/18
NTS-16-030-050 NTS-16-080-050	18.6-16-18.6	20~70*	30~80*	500 500	6/18 6/18
NTS-18-030-050 NTS-18-080-050	20.5-18-20.5	20~70*	30~80*	500 500	8/24 8/24
NTS-20-030-050 NTS-20-080-050	22.7-20-22.7	20~70*	30~80*	500 500	8/24 8/24
NTS-22-030-050 NTS-22-080-050	24.7-22-24.7	20~70*	30~80*	500 500	8/24 8/24
NTS-10-030-090 NTS-10-080-090	14-10-14	20~70*	30~80*	900 900	6/18 6/18
NTS-12-030-090 NTS-12-080-090	16-12-16	20~70*	30~80*	900 900	6/18 6/18
NTS-14-030-090 NTS-14-080-090	16.5-14-16.5	20~70*	30~80*	900 900	6/18 6/18
NTS-16-030-090 NTS-16-080-090	18.6-16-18.6	20~70*	30~80*	900 900	6/18 6/18
NTS-18-030-090 NTS-18-080-090	20.5-18-20.5	20~70*	30~80*	900 900	8/24 8/24
NTS-20-030-090 NTS-20-080-090	22.7-20-22.7	20~70*	30~80*	900 900	8/24 8/24
NTS-22-030-090 NTS-22-080-090	24.7-22-24.7	20~70*	30~80*	900 900	8/24 8/24

\*Increment of 10mm

### HANAROSTENT®

#### Trachea/Bronchium (CCC)



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(mm/Fr)
TCS-10-020-050 TCS-10-080-050	11.5-10-11.5	14~74*	20~80*	500	4/12
TCS-12-020-050 TCS-12-080-050	14.2-12-14.2	13~73*	20~80*	500	4/12
TCS-14-020-050 TCS-14-080-050	15.7-14-15.7	14~74*	20~80*	500	4/12
TCS-16-020-050 TCS-16-080-050	17.5-16-17.5	13~73*	20~80*	500	4/12
TCS-18-020-050 TCS-18-080-050	20-18-20	10~70*	20~80*	500	4/12
TCS-20-020-050 TCS-20-080-050	22-20-22	10~70*	20~80*	500	4/12
TCS-10-020-090 TCS-10-080-090	11.5-10-11.5	14~74*	20~80*	900	4/12
TCS-12-020-090 TCS-12-080-090	14.2-12-14.2	13~73*	20~80*	900	4/12
TCS-14-020-090 TCS-14-080-090	15.7-14-15.7	14~74*	20~80*	900	4/12
TCS-16-020-090 TCS-16-080-090	17.5-16-17.5	13~73*	20~80*	900	4/12
TCS-18-020-090 TCS-18-080-090	20-18-20	10~70*	20~80*	900	4/12
TCS-20-020-090 TCS-20-080-090	22-20-22	10~70*	20~80*	900	4/12

\*Increment of 10mm

### HANAROSTENT®

#### Trachea/Bronchium (NNN)



Model	Stent (mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(mm/Fr)
TNS-12-030-090 TNS-12-080-090	16-12-16	20~70*	30~80*	900	4/12
TNS-14-030-090 TNS-14-080-090	16.5-14-16.5	20~70*	30~80*	900	4/12
TNS-16-030-090 TNS-16-080-090	18.6-16-18.6	20~70*	30~80*	900	4/12

\*Increment of 10mm



※ Sizing and availability varies by country

MEMO

HANAR  STENT®

HANAROSTENT®  
*Broad and Growing Product Selection*

*Türkiye Exclusive Distribütörü*



Güzeltepe Mahallesi Fuar Sokak No: 9/11  
Çankaya / ANKARA

E-Mail: [info@ena-med.com](mailto:info@ena-med.com)  
Tel: 0312 442 95 86

*Karadeniz Bölge Bayi*



Körfez Mahallesi Atatürk Bulvarı No:9  
Atakum / SAMSUN

E - Mail: [diatek@diateksaglik.com.tr](mailto:diatek@diateksaglik.com.tr)  
**Tel: 0362 438 07 50**